

LAMP-Seq Starter Kit

V2.0, 2020/04/16

Jonathan L. Schmid-Burgk^{1,2,3,4,8,13,*}, David Li^{1,2,3,4}, David Feldman⁹, Mikołaj Słabicki^{1,10,11}, Jacob Borrajo^{1,4}, Jonathan Strecker^{1,2,3,4}, Brian Cleary¹, Aviv Regev^{1,5,6,7,12}, Feng Zhang^{1,2,3,4,12}

NOTE: This protocol has not been validated with clinical samples. To facilitate collaborations with interested parties to jointly advance the fight against the current coronavirus pandemic, we have set up a public forum on www.LAMP-Seq.org.

¹ Broad Institute of MIT and Harvard
Cambridge, MA 02142, USA

² McGovern Institute for Brain Research

³ Department of Brain and Cognitive Sciences

⁴ Department of Biological Engineering

⁵ Department of Biology

⁶ Klarman Cell Observatory

⁷ Koch Institute for Integrative Cancer Research

Massachusetts Institute of Technology, Cambridge, MA 02139, USA

⁸ Institute of Clinical Chemistry and Clinical Pharmacology, University Hospital Bonn, 53127
Bonn, Germany

⁹ Department of Biochemistry and Institute for Protein Design, University of Washington,
Seattle, WA 98195, USA

¹⁰ Department of Medical Oncology, Dana-Farber Cancer Institute, Boston, MA 02215, USA

¹¹ Division of Translational Medical Oncology, German Cancer Research Center (DKFZ) and
National Center for Tumor Diseases (NCT), 69120 Heidelberg, Germany

¹² Howard Hughes Medical Institute, Cambridge, MA 02139, USA

¹³ Lead Contact

*Correspondence should be addressed to J.S.-B. (jschmid@broadinstitute.org)

Purpose of this LAMP-Seq Starter Kit

This kit provides detailed instructions for establishing a LAMP-Seq reaction. The purpose of this kit is setting up the method in other labs for further development and optimization.

This protocol has not been validated on clinical samples.

Materials

- 1.5 ml reaction tubes
- PCR strip tubes and caps
- Thermocycler
- Agarose gel equipment
- Ultrapure water
- WarmStart® LAMP Kit (NEB, E1700S)
- NEBNext® High-Fidelity 2X PCR Master Mix (NEB, M0541S)
- 2019-nCoV_N_Positive Control plasmid (IDT, 10006625)
- HiScribe™ T7 High Yield RNA Synthesis Kit (NEB, E2040S)
- Unrelated plasmid (e.g. pX330-U6-Chimeric_BB-CBh-hSpCas9, Addgene #42230)
- Primer C-IVT-fwd (IDT

AAGCTAATACGACTCACTATAGCATACAATGTAACACAAGCTTTTCGG

- Primer C-IVT-rev (IDT)

CTTGATCTTTGAAATTTGGATCTTTGTC

- Primer C-F3 (IDT)
AACACAAGCTTTCGGCAG
- Primer C-B3 (IDT)
GAAATTTGGATCTTTGTCATCC
- Primer C-FIP-AGCTTACT (IDT)
TGCGGCAATGTTTGTAATCAG AGCTTACT CCAAGGAAATTTGGGGAC
- Primer C-BIP (IDT)
CGCATTGGCATGGAAGTCAC TTTGATGGCACCTGTGTAG
- Primer C-LF (IDT)
TTCCTTGTCTGATTAGTTC
- Primer C-LB (IDT)
ACCTTCGGGAACGTGGTT
- Primer PCR-C-fwd-10 (IDT)
ACACTCTTCCCTACACGACGCTCTTCCGATCTAACGCTGAAGCGCTGGGGGC
AAA
- Primer PCR-C-rev-11 (IDT)
TGACTGGAGTTCAGACGTGTGCTCTTCCGATCTGTTTGTAATCAGTTCCTTGTCTG

Protocol

1. Generate template RNA by IVT
 - a. Perform a NEBNext PCR

- i. Primers: C-IVT-fwd and C-IVT-rev
 - ii. Template: 2019-nCoV_N_Positive Control plasmid
 - b. Column-purify and elute in water,
 - c. perform an IVT reaction using HiScribe™ T7 High Yield RNA Synthesis Kit for 12 hours,
 - d. column-purify and elute in water,
 - e. measure concentration using NanoDrop.
2. Make the following plasmid pre-dilutions:
 - a. Make a stock dilution of an unrelated plasmid in water at 5 ng/μl:
 - i. 5 μl plasmid (1 μg/μl) + 995 μl water.
 - b. Create a dilution series:
 - i. Using the plasmid dilution, sequentially dilute template RNA down to four molecules / μl. Use large pipetting volumes if possible (e.g. 2 μl + 998 μl).
3. Prepare a 10x LAMP primer mix:
 - a. 16 μl C-FIP-AGCTTACACT (100 μM),
 - b. 16 μl C-BIP (100 μM),
 - c. 2 μl C-F3 (100 μM),
 - d. 2 μl C-B3 (100 μM),
 - e. 4 μl C-LF (100 μM),
 - f. 4 μl C-LB (100 μM),
 - g. 56 μl water,
 - h. mix by pipetting.
4. Perform barcoded LAMP reactions:

- a. Prepare a reaction mix:
 - i. 312.5 μ l 2x LAMP Master Mix (NEB),
 - ii. 62.5 μ l 10x LAMP primer mix,
 - iii. 125 μ l water,
 - iv. 125 μ l template RNA (four molecules / μ l RNA in 5 ng/ μ l plasmid DNA),
 - v. mix by pipetting.
 - b. Aliquot 25 μ l reaction mix into each of 24 tubes in PCR strips,
 - c. incubate the LAMP reactions in a thermocycler for 40 minutes at 65 °C, and subsequently for 10 minutes at 95 °C.
5. Dilute the LAMP reactions 100-fold in PCR strips:
- a. 2 μ l LAMP reaction + 198 μ l water.
6. Perform PCR reactions:
- a. Prepare a reaction mix:
 - i. 150 μ l NEBNext 2x Master Mix (NEB),
 - ii. 1.5 μ l PCR-C-fwd-10 (100 μ M),
 - iii. 1.5 μ l PCR-C-rev-11 (100 μ M),
 - iv. 147 μ l water,
 - v. mix by pipetting.
 - b. Aliquot 9 μ l reaction mix into each of 24 tubes in PCR strips,
 - c. to each reaction, add 1 μ l of diluted LAMP reaction,
 - d. mix by pipetting,
 - e. run the PCR reaction in a thermocycler with the following protocol:
 - i. 98 °C 3 minutes,

- ii. 98 °C 20 seconds (15 cycles),
- iii. 65 °C 20 seconds (15 cycles),
- iv. 72 °C 30 seconds (15 cycles),
- v. 72 °C 3 minutes,
- vi. 4 °C hold.

7. Run 5 μ l of each PCR reaction on a 1% agarose gel.

Expected outcome

Positive reactions are expected to yield a single strong band at an estimated size of 184 bp. At 20 molecules per reaction, around 10 out of 24 reactions are expected to be detected as positive.

